HÆGER

Body Analyser scale with app / Balança de análise corporal com app
Báscula analizador corporal con app / Balance d'analyse corporelle avec
application / Ζυγαριά αναλυτή σώματος με εφαρμογή

Download
AiLink app

InteliScale * BS-DIG.014A

Download
Multi-language file

CE RoHS

User instructions / Instruções de uso
Instrucciones de uso / Mode d'emploi
Οδηγίες χρήστη

- Keep the scale out of water.
- Mantenha a balança afastada de água.
- Mantenga la balanza alejada del agua.
- Gardez la balance loin de l'eau.
- Κρατήστε την κλίμακα έξω από το νερό.

- Keep the scale with care.
- Transporte a balança com cuidado.
- Transporte la balanza cuidadosamente.
- Transportez la balance soigneusement.
- Κρατήστε την κλίμακα με προσοχή.

- Don't put it on uneven floor while weighing.
- Não coloque sobre piso irregular durante a pesagem.
- No lo ponga en el suelo irregular, mientras que está pesando.
- Ne le placez pas sur un sol irrégulier tout en pesant.
- Μην το βάλετε σε άνισο πάτωμα ενώ ζυγίζει.

- Don't step on the edge of the scale.
- Não pise nas extremidades da balança.
- No pise en el borde de la escala.
- Ne marchez pas sur les bords de l'échelle.
- Μην βγαίνουμε στα άκρα της κλίμακας.

English

Dear customer

Thank you for having chosen a HÆGER product. The HÆGER products have been produced to think about the welfare of the consumer privileging the most raised standards of quality, functionality and design. We are sure you will be happy with this appliance.

We assume that the user is familiar with the common procedures of handling household appliances.

To ensure the scale continues to deliver peak performance, please read these operation instructions carefully and thoroughly prior to initial usage.

General Information

Your scale is a precision electronic instrument. Treat your scale with care. The scale is durable, but misuse or abuse may render the electronic sensors inoperable.

1. Read all instructions before operating the scale.
2. Follow the basic safety precautions for all office equipment when using this scale.
3. The manufacturer does not warrant this unit to be waterproof. Do not expose this product to rain or moisture.
4. Do not use the scale under sun or keep it away from heating appliances such as oven, stove, etc.
5. This scale does not contain any user-serviceable parts. Do not disassemble.
6. Use the scale for its intended purpose only.

Note: *The display value may be affected by electromagnetic influences, e.g., when a radio is operated in the immediate vicinity of the device. The product can be used for its intended purpose again when the interference is discontinued but may need*

to be switched on again to re-start.

! WARNING

- DO NOT use with medical implants such as pacemakers.
- DO NOT stand on the edge of the scale or jump on it.
- DO NOT overload the scale (MAX. 180kg / 396lb).
- DO NOT drop scale or droop any objects on it as this may damage the sensors.
- DO NOT immerse the scale in water or use chemical cleaning agents. Clean the scale with slightly damp cloth.
- NOT recommended for pregnant women.
- NOT recommended for infants, toddlers, and children under 10 years of age.
- Any measurement obtained using this device is for reference only and should not be considered as a medical opinion.
- Consult your primary care doctor or physician before making changes to your diet, exercise plan or physical activities.
- Always place the scale on a hard, dry and flat surface before measurement.
- Make sure your feet are dry before stepping on the scale.
- Connect the batteries in accordance with the correct polarities indicated.
- Remove batteries if the scale is not be used for a prolonged period of time.
- The user under 6 years old, only weight function, no BMI and body fat, etc. data.
- The user under 10 years old, only weight & BMI function, no body fat, etc. data.

Features

- Maximum Weight: 180kg/396lb
- Up to 8 different users in the app
- LED screen
- Unit: kg/lb
- Battery: 3xAAA batteries — Not included
- Frequency Range: 2402~2480MH
- Maximum transmitted radio frequency

- power: -0.34dBm
- Weight tolerance: 50kg±300g; 100kg±400g; 150kg±500g
 - Automatically turns off
 - Low battery indicator
 - Wireless technology
 - App compatible with App Store and Google Play (Trademarks)
 - Measurement of Body Mass Index
 - Measurement of bone mass and muscle mass
 - Visceral fat index
 - Body fat rate
 - Other measurements (consult the application)

Battery Installation

1. Turn the entire scale upside down and reset it securely on a firm, flat surface.
2. Remove the battery door.
3. Install 3x1.5V AAA battery in the battery compartment.
4. Close the battery door.

Battery Status

- When the low battery symbol **LD** appears on the LCD screen, it indicates that the battery is low. Open the battery compartment and replace the old batteries with new ones.
- **Used batteries must not be disposed of as household waste. Do not incinerate batteries because they might explode at high temperature. Contact your local authority for information concerning reclamation and disposal of batteries.**
- **When replacing batteries, always use brand new battery of the same type.**

Download Application

1. Download the application via QRCode in the cover page. (The application is compatible with Android and IOS system)

2. Register the user in the application.
3. Update profile and requested data.
4. Turn on Bluetooth on the mobile phone.
5. Open the application and select the device select "Body Fat Scale (Bluetooth).
6. Pair the mobile phone with the scale.
7. Position yourself on the scale and wait for the data to synchronize.

Operation

1. Place the scale on a flat, hard surface before weighing.
2. Keep the surface of the scale and feet dry during weighing. Slippery surface when wet.
3. The scale switches on and off automatically.
4. Download the application using the QRCode found on the packaging, in the manual and on the scale itself.
5. Turn on the mobile phone's Bluetooth to connect the application with the scale. If you don't turn on Bluetooth, just set on top of the scale and it will only show the weight.
6. Register and follow all the steps requested in the application.
7. Replace the batteries immediately when the low battery indicator appears.
8. Also check that the battery polarity is correct.
9. When **Err** appears on the screen, it means that the scale is overloaded.

Usage tips

- The scales are designed for a maximum weight of 180kg (396lb).
- Weigh yourself without any clothes or shoes on before meals and always at the same hour of the day. The best time to

weigh yourself is when you wake up in the morning.

- The balance rounds off weight to the closest 100 g.
- If you weigh yourself twice and two different values are given, it means that your weight is comprised between both values.

Care and Maintenance

- Do not leave any objects on top of the scale when not in use.
- Keep the surface of the balance and feet dry during weighing. Slippery surface when wet.
- Make sure you are in the centre of the scale to avoid tipping over.
- Avoid dropping or hitting the scale. It will damage it.
- Wipe the surface of the balance with a clean, dry cloth.
- Do not wash the balance with water, the balance mechanism will be damaged if submerged in water.
- If the scale will not be used for long periods, it is recommended to remove the batteries to avoid possible damage to the balance.
- Waste electrical products must not be disposed of with household waste.
- Contact your local authority or dealer for recycling advice.

Technical Specifications

Model:	BS-DIG.014A
Maximum Capacity:	180kg (396lb)
Graduation:	100g/0.2 lb
Power Source:	3x1.5V AAA batteries

Conformity

This device has been tested according to all relevant current CE guidelines, such as

electromagnetic compatibility Directive 2014/30/EU and Radio Directive RED 2014/53/EU, and has been constructed in accordance with the latest safety regulations.

Disposal – Environment policy

Disposal

This appliance is marked according to the European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE).

By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

This symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead, it shall be

handed over to the applicable collection point for the recycling of electrical and electronic equipment.

Disposal must be carried out in accordance with local environmental regulations for waste disposal. For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or the shop where you purchased the product.

Disposal of used batteries

Batteries contain environmentally damaging compounds and, therefore, do not belong in standard communal waste. Deliver used batteries for proper disposal at locations intended for their collection.

Guarantee

Warranty Terms and Conditions

1. The product is guaranteed for the period stipulated by the legislation where it is placed on the market, starting from the date of purchase, against manufacturing defects and/or workmanship. The lack of conformity manifested in this period is presumed to exist at

- the date of purchase. After this period, it is for the customer to prove that the lack of conformity already existed at the time of delivery.
2. If, during the warranty period, the product is defective in normal use, the buyer should contact the selling establishment and send the appliance to the place indicated by the seller.
 3. The warranty is only valid upon presentation of the invoice proving the purchase and the completed warranty certificate (indicating the date of purchase, the name of the dealer, the reference of the model, and it is also recommended to indicate the serial number and lot number).
 4. The importer/seller reserves the right to refuse warranty assistance in the event that such information has been deleted or altered after the original purchase of the product.
 5. The responsibility of the importer/seller includes in particular the costs of repairing and/or replacing the unit covered by the warranty with the reservation of the right to replace it with an equivalent product, in cases where it is not possible to repair it.
 6. The warranty does not apply to problems that are not directly related to defects in material, design or workmanship.
 7. This product is an appliance and is intended for home use only. The warranty will be void if the product is used for professional purposes.
 8. The warranty does not cover damage caused by misuse, improper installation, spillage, external factors, or intentional damage.
 9. The warranty does not cover damage resulting from falls, use of excessive force, blows, exposure to extreme environmental conditions or deterioration caused by the normal use of plastic parts or keyboards, as well as the use of batteries other than those specified in this manual.
 10. The warranty does not cover installation and programming by the dealer. In particular, no claim for installation and free programming by the reseller may be accepted.
 11. The warranty does not cover deficiencies that insignificantly impair the operation of the appliance. Other costs, such as installation, transportation and travel of the technician are expressly excluded from the warranty.
 12. The applicability of the warranty presupposes that the appliance is sent to the seller or place indicated by him, complete and in well-protected packaging (if possible, in its original packaging) and accompanied by its sheet or warranty card duly filled out and proof of purchase.
 13. If after verification, it is found that there are no reasons for the claim, or that the product is not defective, the inherent costs will be imputed to the customer, and the dealer is authorized to charge those costs to the customer.
 14. The warranty will be void when it shows signs that an unauthorized person has attempted to carry out repairs, modifications or replacement of parts on the appliance.
 15. Repairs carried out after the end of the warranty period are subject to costs.
 16. This warranty does not affect your legal rights that you may have as a consumer under applicable national law governing the purchase of products where the product is sold.
 17. **Warranty Period:**
 - European Union countries** - European Directives (EU) 2019/771 and (EU) 2019/770 apply on certain aspects of contracts for the purchase and sale of movable property and transposed into the national laws of their countries.
 - Other countries** - national legislation on consumer rights in their country applies. In the case of no official legislation, the guarantee shall be applied at the discretion of the importer placing the product on the market or the seller establishment.

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HAEGERTEC, S.A.